


SCAKIT

Sinus Crestal Approach



Drill Sinus Cortical Wall with No Fear! Crestal Approach without Membrane Damage

SCA kit

SCA kit is an effective tool for crestal sinus elevation by perforating inferior cortical wall and lifting sinus membrane without any rupture or membrane tearing. (1,200 RPM with irrigation)

SCA kit's 3S

SAFETY

S-reamer's head with a special blade structure is designed to perforate the sinus inferior cortical bone without sinus membrane tearing and it makes thin bone disk between membrane and reamer. Since the reamer does not touch sinus membrane directly, so it is safer.

SPEEDY

With 1,200 RPM high speed drilling with enough irrigation perforates sinus cortical wall safely and rapidly without any sinus membrane damage.

SIMPLE

SCA kit is designed user-friendly to minimize the operators' experience gap and sinus lifting techniques, therefore all users can conduct sinus lift easily.

