


Neo Sinus All kit




Neo Sinus All kit



Neo Biotech
 Satisfaction to Dentists

Composition of Neo Sinus All kit


Crestal Approach



Lateral Approach

Neo Biotech
Satisfaction to Dentists

Composition of Neo Sinus All kit



Neo Biotech
Satisfaction to Dentists

Crestal Approach Part




Comparison – Crestal part

SCA Kit



Ø2.0

Sinus All Kit



Ø2.2

Ø2.8

Ø3.2

Ø3.7

Ø4.0

Can complete drilling for the implantation with Sinus All Kit

Comparison – Crestal part

SCA Kit




Ø2.4

Ø2.8

Ø3.2

Ø3.6

Sinus All Kit



Ø2.8

Ø3.2

Ø3.6

Ø3.9

S-Reamer size is changed

Comparison – Crestal part

SCA Kit



Sinus All Kit



- ✓ 1mm, 12mm stopper
- ✓ Connection type
- ✓ Irrigation
- ✓ Colorful

Comparison – Depth Gauge

SCA Kit




Sinus All Kit



2 diameters

Comparison – Rotary elevator


SCA Kit



Inserter

spreader

Sinus All Kit



Ø2.8


Ø3.2

Ø3.6


Ø3.9

Rotary elevator

Comparison – Rotary elevator




Sinus All Kit



Rotary elevator

Comparison – Bone condensor

SCA Kit



2mm



3mm



4mm



5mm



6mm



7mm



8mm



9mm




10mm



11mm

Sinus All Kit



1mm



2mm



3mm



4mm



5mm



6mm



7mm



8mm



9mm



10mm



11mm




12mm

Addition – Drill handle



- ✓ Hand motor
- ✓ Can be mounted contra angle connection
- ✓ Stable operation

Addition – driver, ratchet




Fixture driver



Torque Ratchet

Addition – Aqua system



- ✓ Sinus membrane detach with water
- ✓ 3cc syringe
- ✓ Stable lifting (“Gauge” control the pressure)

Addition – Aqua system



Lateral Approach Part



Comparison – LS reamer



SLA Kit



Sinus All Kit



- ✓ Simple
- ✓ Improve cutting ability
- ✓ Stopper

Comparison – C guide



SLA Kit




Sinus All Kit




Comparison – Elevator



SLA Kit



Sinus All Kit



Addition – Side cutter



- ✓ Widening window
- ✓ Stopper

Addition – Side cutter



- ✓ Widening window
- ✓ Stopper